


Deltalab Natural Polyethylene Caps for 16mm Conical Tubes (with Flange Plug)

Product Code: 305502

Brand: Deltalab

Ideal for capping and recapping 16mm conical tubes, Deltalab's Natural Polyethylene Caps for 16mm Conical Tubes (with Flange Plug) present significant advantages over other flange caps on the market, which are:

- A larger size, providing a larger surface for finger rest
- Flanges have a large size than on other caps. More space between each flange increases the cap's sealing qualities
- A small internal cavity designed to prevent possible drop spillage
- A raised top for easy-on, easy-off use and to avoid cap deformation from repeated removal
- A ribbed surface.



Specifications:

Colour	Natural
Internal Diameter (mm)	14.50
External Diameter (mm)	17.05
Total Length (mm)	21.77
Material	Polyethylene
Sterility	Non-Sterile
Quantity	2500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com