

Deltalab Natural Polyethylene Caps for 12mm Tubes (Ribbed)

Product Code: 305802

Brand: Deltalab

Ideal for capping and recapping 12mm blood collection tubes, Deltalab's Natural Polyethylene Caps for 12mm Tubes (Ribbed) are designed to be placed and removed securely. Each cap has a small internal cavity designed to prevent possible drop spillage.

Ergonomical, these leakproof polyethylene caps are ribbed for better handling while using gloves. The top of the caps are raised for easy-on, easy-off use and to avoid cap deformation from repeated removal.

Specifications:

Colour	Natural
Internal Diameter (mm)	10.60
External Diameter (mm)	15
Total Length (mm)	15.30
Material	Polyethylene
Sterility	Non-Sterile
Quantity	2000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com