

## Deltalab Natural Polyethylene Caps for 12mm Tubes (Ribbed)

**Product Code:** 305802

**Brand:** Deltalab

Ideal for capping and recapping 12mm blood collection tubes, Deltalab's Natural Polyethylene Caps for 12mm Tubes (Ribbed) are designed to be placed and removed securely. Each cap has a small internal cavity designed to prevent possible drop spillage.

Ergonomical, these leakproof polyethylene caps are ribbed for better handling while using gloves. The top of the caps are raised for easy-on, easy-off use and to avoid cap deformation from repeated removal.



### Specifications:

<b>Colour</b>	Natural
<b>Internal Diameter (mm)</b>	10.60
<b>External Diameter (mm)</b>	15
<b>Total Length (mm)</b>	15.30
<b>Material</b>	Polyethylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	2000

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)