

Deltalab Green Polyethylene Caps for 12mm Tubes (Ribbed)

Product Code: 305807

Brand: Deltalab

Ideal for capping and recapping 12mm blood collection tubes, Deltalab's Green Polyethylene Caps for 12mm Tubes (Ribbed) are designed to be placed and removed securely. Each cap has a small internal cavity designed to prevent possible drop spillage.

Ergonomical, these leakproof polyethylene caps are ribbed for better handling while using gloves. The top of the caps are raised for easy-on, easy-off use and to avoid cap deformation from repeated removal.



Specifications:

| | |
|-------------------------------|--------------|
| Colour | Green |
| Internal Diameter (mm) | 10.60 |
| External Diameter (mm) | 15 |
| Total Length (mm) | 15.30 |
| Material | Polyethylene |
| Sterility | Non-Sterile |
| Quantity | 2000 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com