

Deltalab Red Polyethylene Caps for 16mm Tubes (Ribbed)


Product Code: 307006

Brand: Deltalab

Ideal for capping and recapping tubes of different sizes, Deltalab's Red Polyethylene Caps for 16mm Tubes (Ribbed) are designed to be placed and removed securely. Each cap is ribbed for better handling while using gloves.

They have four internal rims that can be clamped firmly on the tube. These reversible caps can be used as a watertight seal or inverted in its recap position as a simple closure when no leakproof seal is required.

Please note that these caps can be used on Deltalab's 17ml Round Bottom Polystyrene Tubes (16 x 150mm).



Specifications:

Colour	Red
Internal Diameter (mm)	14.80
External Diameter (mm)	19.20
Total Length (mm)	25.35
Material	Polyethylene
Sterility	Non-Sterile
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com