


Deltalab Natural Polyethylene Caps for 12 - 13mm Tubes (with Flange Plug)

Product Code: 307102

Brand: Deltalab

Ideal for capping and recapping 12 - 13mm conical tubes, Deltalab's Natural Polyethylene Caps for 12 - 13mm Tubes (with Flange Plug) present significant advantages over other flange caps on the market, which are:

- A larger size, providing a larger surface for finger rest
- Flanges have a large size than on other caps. More space between each flange increases the cap's sealing qualities
- A small internal cavity designed to prevent possible drop spillage
- A raised top for easy-on, easy-off use and to avoid cap deformation from repeated removal
- A ribbed surface.



Specifications:

Colour	Natural
Internal Diameter (mm)	11.25
External Diameter (mm)	13.30
Total Length (mm)	17.40
Material	Polyethylene
Sterility	Non-Sterile
Quantity	2000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com