


## Deltalab Polyethylene Caps for 12 - 13mm Tubes (with Flange Plug)

**Product Code:** 307108

**Brand:** Deltalab

Ideal for capping and recapping 12 - 13mm conical tubes, Deltalab's Mauve Polyethylene Caps for 12 - 13mm Tubes (with Flange Plug) present significant advantages over other flange caps on the market, which are:

- A larger size, providing a larger surface for finger rest
- Flanges have a large size than on other caps. More space between each flange increases the cap's sealing qualities
- A small internal cavity designed to prevent possible drop spillage
- A raised top for easy-on, easy-off use and to avoid cap deformation from repeated removal
- A ribbed surface.



### Specifications:

|                               |              |
|-------------------------------|--------------|
| <b>Colour</b>                 | Mauve        |
| <b>Internal Diameter (mm)</b> | 11.25        |
| <b>External Diameter (mm)</b> | 13.30        |
| <b>Total Length (mm)</b>      | 17.40        |
| <b>Material</b>               | Polyethylene |
| <b>Sterility</b>              | Non-Sterile  |
| <b>Quantity</b>               | 2000         |

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)