


## Deltalab Polyethylene Caps for 15 - 16mm Tubes (with Flange Plug)

**Product Code:** 307205

**Brand:** Deltalab

Ideal for capping and recapping 15 - 16mm conical tubes, Deltalab's Blue Polyethylene Caps for 15 - 16mm Tubes (with Flange Plug) present significant advantages over other flange caps on the market, which are:

- A larger size, providing a larger surface for finger rest
- Flanges have a large size than on other caps. More space between each flange increases the cap's sealing qualities
- A small internal cavity designed to prevent possible drop spillage
- A raised top for easy-on, easy-off use and to avoid cap deformation from repeated removal
- A ribbed surface.



### Specifications:

<b>Colour</b>	Blue
<b>Internal Diameter (mm)</b>	15.40
<b>External Diameter (mm)</b>	17
<b>Total Length (mm)</b>	21.80
<b>Material</b>	Polyethylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	2500

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)