

Deltalab Natural Polyethylene Caps for 16mm Conical Tubes (Ribbed)

Product Code: 308202

Brand: Deltalab

Offering one of the largest chambers on the market, Deltalab's Natural Polyethylene Caps for 16mm Conical Tubes (Ribbed) are designed to be placed and removed securely on 16mm conical tubes. With a leakproof closing, these caps allow an almost complete filling of the tube and eliminates the need for a fill-in line.

They are ribbed for better handling while using gloves.

The inner part of the caps, long-sized, consist of three rings to ensure a tight, safe leakproof closing. The top of the caps are raised for easy-on, easy-off use and to avoid cap deformation from repeated removal. They are ideal for use with urine.

To avoid possible liquid spillage on the bench top, leave the caps to stand with flanges up.

Specifications:

Colour	Natural
Internal Diameter (mm)	14.55
External Diameter (mm)	20.45
Total Length (mm)	18.50
Material	Polyethylene
Sterility	Non-Sterile
Quantity	2500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

