

Medline Polystyrene Clear Dipper Samplers with Handles (30ml)

Product Code: 309277

Brand: Medline Scientific

Medline's Polystyrene Clear Dipper Samplers with Handles (30ml) allows laboratories to carry out simple sample collections and subsequent transportation. The handle neatly snaps off after sample collection from the vat for the ease of transportation.

Each dipper is supplied with a flow-seal screw cap to ensure leak-proof closure.

Specifications

Colour	Clear
Volume (ml)	30
Handle Length (mm)	385
Material	Polystyrene
Inner Pack	1
Quantity	50

