

## Medline Polystyrene Containers with Labels (60ml)

**Product Code:** 309501P


**Brand:** Medline Scientific

Manufactured from medical grade virgin polystyrene, Medline's Polystyrene Containers with Labels (60ml) have a white screw cap that allows leaks to be avoided.

All containers are CE marked in accordance with the In Vitro Device Directive.

### Specifications:

<b>Capacity (ml)</b>	60
<b>Material</b>	Polystyrene
<b>Quantity</b>	300



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)