

## Medline Polystyrene Universal Containers with Spoons & Printed Labels (30ml)

**Product Code:** 309602E

**Brand:** Medline Scientific

Made of polystyrene and supplied sterile, Medline's Polystyrene Universal Containers with Spoons & Printed Labels (30ml) are suitable for collecting faecal or semi-solid samples through the use of their spoons, which are sealed into their screw caps.

### Specifications:

<b>Capacity (ml)</b>	30
<b>Material(s)</b>	Polystyrene
<b>Cap Colour</b>	White
<b>Sterility</b>	Sterile
<b>Quantity</b>	400



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)