

Medline Polystyrene Blue Dipper Samplers with Handles (100ml)

Product Code: 309688

Brand: Medline Scientific

Medline's Polystyrene Blue Dipper Samplers with Handles (100ml) allow laboratories to carry out simple sample collections and subsequent transportation. The handle neatly snaps off after sample collection from the vat for the ease of transportation.

Each dipper is supplied with a flow-seal screw cap to ensure leak-proof closure.

Specifications

Colour	Blue
Volume (ml)	100
Handle Length (mm)	330
Material	Polystyrene
Inner Pack	1
Quantity	50



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com