

## Ceti Triton II Multi-Observation Trinocular Compound Microscope with Two Teaching Heads

**Product Code:** 3100.9100M

**Brand:** Ceti

The Ceti Triton II Multi-Observation Trinocular Compound Microscope with Two Teaching Heads is a research microscope with 100W illumination and an two additional teaching heads.



### Specifications

<b>Optical Head</b>	Trinocular Ergonomic 0 - 45°
<b>Additional Number of Optical Heads</b>	2
<b>Eye Piece</b>	10x/20mm wide-field (pair)
<b>Objectives</b>	4x, 10x, 40x, 100x E-plan Plan Achromatic Infinity Corrected
<b>Condenser</b>	Swing-out condensing lens - Brightfield NA 1.25
<b>Stage</b>	Size 185 x 142mm with built-in mechanical stage and slide clip
<b>Focusing</b>	Coaxial course and fine
<b>Köhler Illumination (variable)</b>	Halogen 100W

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)