

Ceti Triton II Multi-Observation Trinocular Compound Microscope with Three Teaching Heads

Product Code: 3100.9200M

Brand: Ceti

The Ceti Triton II Multi-Observation Trinocular Compound Microscope with Three Teaching Heads is a research microscope with 100W illumination and three additional teaching heads.



Specifications

Optical Head	Trinocular Ergonomic 0 - 45°
Additional Number of Optical Heads	3
Eye Piece	10x/20mm wide-field (pair)
Objectives	4x, 10x, 40x, 100x E-plan Plan Achromatic Infinity Corrected
Condenser	Swing-out condensing lens - Brightfield NA 1.25
Stage	Size 185 x 142mm with built-in mechanical stage and slide clip
Focusing	Coaxial coarse and fine
Köhler Illumination (variable)	Halogen 100W

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com