


## Deltalab 48ml Pipetting Reservoirs (8 channel)

**Product Code:** 330089

**Brand:** Deltalab

Made of polystyrene, Deltalab's 48ml Pipetting Reservoirs (8 channel) provide a strong clearance for various types of pipettors. They are perfect for carrying out long multi-channel pipetting into microwell plates. Each reservoir is sterile by radiation.



## Specifications

<b>Dimensions (mm)</b>	83 x 126 x 13
<b>Description</b>	8 channel, 48ml
<b>Material</b>	Polystyrene
<b>Sterility</b>	Sterile
<b>Quantity</b>	10

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)