

## Deltalab Medical Grade Silicone Tubes (0.5 x 1mm)

**Product Code:** 350051

**Brand:** Deltalab

Deltalab's Medical Grade Silicone Tubes (0.5 x 1mm) can be used for a variety of medical, pharmaceutical, industrial, and clinical applications. These translucent tubes are highly mechanical and chemically resistant. They have a chemical composition similar to quartz and glass and share their excellent properties.

Each tube is autoclavable and made of materials suitable for alimentary use.



## Specifications

<b>Specific Gravity</b>	1.2 ± 0.1
<b>Toughness</b>	40 - 80 Sh
<b>Resistance to Traction</b>	77 - 95kg/cm <sup>2</sup>
<b>Extension</b>	275 - 780%
<b>Resistance to Tearing</b>	20 - 50kg/cm <sup>2</sup>
<b>Temperature Range</b>	-50°C to 200°C (beaks of 250°C)
<b>Dimensions (Ø int. x Ø ext.)</b>	0.5 x 1mm
<b>Quantity</b>	100

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)