

Deltalab Medical Grade Silicone Tubes (0.5 x 1mm)

Product Code: 350051

Brand: Deltalab

Deltalab's Medical Grade Silicone Tubes (0.5 x 1mm) can be used for a variety of medical, pharmaceutical, industrial, and clinical applications. These translucent tubes are highly mechanical and chemically resistant. They have a chemical composition similar to quartz and glass and share their excellent properties.

Each tube is autoclavable and made of materials suitable for alimentary use.

Specifications

Specific Gravity	1.2 + 0.1
Toughness	40 - 80 Sh
Resistance to Traction	77 - 95kg/cm ²
Extension	275 - 780%
Resistance to Tearing	20 - 50kg/cm ²
Temperature Range	-50°C to 200°C (beaks of 250°C)
Dimensions (Ø int. x Ø ext.)	0.5 x 1mm
Quantity	100

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com