

## Dentalab Medical Grade Silicone Tubes (5 x 7mm)

**Product Code:** 350057

**Brand:** Dentalab

Dentalab's Medical Grade Silicone Tubes (5 x 7mm) can be used for a variety of medical, pharmaceutical, industrial, and clinical applications. These translucent tubes are highly mechanical and chemically resistant. They have a chemical composition similar to quartz and glass and share their excellent properties.

Each tube is autoclavable and made of materials suitable for alimentary use.

### Specifications

|                                     |                                 |
|-------------------------------------|---------------------------------|
| <b>Specific Gravity</b>             | 1.2 + 0.1                       |
| <b>Toughness</b>                    | 40 - 80 Sh                      |
| <b>Resistance to Traction</b>       | 77 - 95kg/cm <sup>2</sup>       |
| <b>Extension</b>                    | 275 - 780%                      |
| <b>Resistance to Tearing</b>        | 20 - 50kg/cm <sup>2</sup>       |
| <b>Temperature Range</b>            | -50°C to 200°C (beaks of 250°C) |
| <b>Dimensions (Ø int. x Ø ext.)</b> | 5 x 7mm                         |
| <b>Quantity</b>                     | 15                              |

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)