

Dentalab Medical Grade Silicone Tubes (10 x 14mm)

Product Code: 351014

Brand: Dentalab

Dentalab's Medical Grade Silicone Tubes (10 x 14mm) can be used for a variety of medical, pharmaceutical, industrial, and clinical applications. These translucent tubes are highly mechanical and chemically resistant. They have a chemical composition similar to quartz and glass and share their excellent properties.

Each tube is autoclavable and made of materials suitable for alimentary use.



Specifications

| | |
|-------------------------------------|---------------------------------|
| Specific Gravity | 1.2 ± 0.1 |
| Toughness | 40 - 80 Sh |
| Resistance to Traction | 77 - 95kg/cm ² |
| Extension | 275 - 780% |
| Resistance to Tearing | 20 - 50kg/cm ² |
| Temperature Range | -50°C to 200°C (beaks of 250°C) |
| Dimensions (Ø int. x Ø ext.) | 10 x 14mm |
| Quantity | 15 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com