

Deltalab Ultra Clear Polypropylene Tubes (12 x 88mm)

Product Code: 400400

Brand: Deltalab


Manufactured from high technology, Deltalab's Ultra Clear Polypropylene Tubes (12 x 88mm) feature new generation moulds to ensure reproducibility. The tubes are made in one piece to ensure uniformity and dimensional accuracy.

They are autoclavable up to 121°C.

Each tube comes with graduation rings graduated at 1, 2.5, and 5ml.

Specifications:

Dimensions (mm)	12 x 88
Volume (ml)	6
Internal Diameter (mm)	10.3
External Diameter (mm)	11.9
Total Length (mm)	86.9
Maximum Volume (ml)	5.8
Material	Polypropylene
Sterility	Non-Sterile
Quantity	4000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com