

Deltalab Ultra Clear Polypropylene Tubes (11 x 55mm)

Product Code: 400500


Brand: Deltalab

Manufactured from high technology, Deltalab's Ultra Clear Polypropylene Tubes (11 x 55mm) feature new generation moulds to ensure reproducibility. The tubes are made in one piece to ensure uniformity and dimensional accuracy.

They are autoclavable up to 121°C.

Specifications:

Dimensions (mm)	11 x 55
Volume (ml)	3
Internal Diameter (mm)	9.9
External Diameter (mm)	11.4
Total Length (mm)	53.9
Maximum Volume (ml)	3.4
Material	Polypropylene
Sterility	Non-Sterile
Quantity	4000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com