

Deltalab Ultra Clear Polypropylene Tubes (13 x 75mm)

Product Code: 400700

Brand: Deltalab

Manufactured from high technology, Deltalab's Ultra Clear Polypropylene Tubes (13 x 75mm) feature new generation moulds to ensure reproducibility. The tubes are made in one piece to ensure uniformity and dimensional accuracy.

They are autoclavable up to 121°C.

Specifications:

Dimensions (mm)	13 x 75
Volume (ml)	5
Internal Diameter (mm)	10.9
External Diameter (mm)	12.7
Total Length (mm)	73.7
Maximum Volume (ml)	6.4
Material	Polypropylene
Sterility	Non-Sterile
Quantity	4000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com