

## Deltalab Ultra Clear Polypropylene Tubes (12 x 75mm)

**Product Code:** 400800

**Brand:** Deltalab


Manufactured from high technology, Deltalab's Ultra Clear Polypropylene Tubes (12 x 75mm) feature new generation moulds to ensure reproducibility. The tubes are made in one piece to ensure uniformity and dimensional accuracy.

They are autoclavable up to 121°C.

Each tube has graduation rings graduated at 1, 2.5, and 4ml.

### Specifications:

<b>Dimensions (mm)</b>	12 x 75
<b>Volume (ml)</b>	5
<b>Internal Diameter (mm)</b>	10.3
<b>External Diameter (mm)</b>	11.7
<b>Total Length (mm)</b>	73.9
<b>Maximum Volume (ml)</b>	5.3
<b>Material</b>	Polypropylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	4000



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)