

Deltalab Ultra Clear Polypropylene Tubes (16 x 95mm)

Product Code: 400900

Brand: Deltalab


Manufactured from high technology, Deltalab's Ultra Clear Polypropylene Tubes (16 x 95mm) feature new generation moulds to ensure reproducibility. The tubes are made in one piece to ensure uniformity and dimensional accuracy.

They are autoclavable up to 121°C.

Each tube has graduation rings graduated with 2.5, 5, and 10ml.

Specifications:

Dimensions (mm)	16 x 95
Volume (ml)	10
Internal Diameter (mm)	13.7
External Diameter (mm)	15.5
Total Length (mm)	93.7
Maximum Volume (ml)	11.9
Material	Polypropylene
Sterility	Non-Sterile
Quantity	2500



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com