

Deltalab 12ml Conical Polypropylene Tubes (16 x 102mm)

Product Code: 401200


Brand: Deltalab

Deltalab's 12ml Conical Polypropylene Tubes (16 x 102mm) can be used for cell centrifugation and other cell culture applications.

Each conical tube have internal graduations and an external rim.

Specifications:

Autoclavable	Yes
Dimensions (mm)	16 x 102
Material	Polypropylene
Volume (ml)	12
Internal Diameter (mm)	14.4
External Diameter (mm)	18.8/16.2
Total Length (mm)	99.8
Max. Volume (ml)	12.5
Quantity	2500



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com