

## Deltalab EUROTUBO 12ml Round Bottom Screw Cap Tubes (Polypropylene - Individually Wrapped)

**Product Code:** 301402

**Brand:** Deltalab

Made of transparent polypropylene, Deltalab's EUROTUBO 12ml Round Bottom Screw Tubes (Polypropylene - Individually Wrapped) feature a green cap made of high density polyethylene. The external skirt allows the tubes to remain free-standing. Each tube is supplied screwed.

### Specifications:

Description	Polypropylene - Individually Wrapped
Recommended Volume (ml)	12
Dimensions (mm)	15 x 102
Autoclavable	No
Sterility	Sterile (ethylene oxide)
Max. Volume (ml)	14.2
External Cap (Ø mm)	20.9
External Tube (Ø mm)	16.2
Internal Tube (Ø mm)	14.3
Length with Cap (Ø mm)	102.9
Length without Cap (Ø mm)	101.5
Quantity	1500

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)