


Deltalab EUROTUBO 12ml Round Bottom Screw Tap Tubes (High Transparency Polypropylene)

Product Code: 401403

Brand: Deltalab

Made of transparent polypropylene, Deltalab's EUROTUBO 12ml Round Bottom Screw Tubes (High Transparency Polypropylene) feature a green cap made of high density polyethylene. The external skirt allows the tubes to remain free-standing. Each tube is supplied screwed.



Specifications:

Description	High Transparency Polypropylene
Recommended Volume (ml)	12
Dimensions (mm)	15 x 102
Autoclavable	Yes
Sterility	Non-Sterile
Max. Volume (ml)	14.2
External Cap (Ø mm)	20.9
External Tube (Ø mm)	16.2
Internal Tube (Ø mm)	14.3
Length with Cap (Ø mm)	102.9
Length without Cap (Ø mm)	101.5
Quantity	1500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com