

## Deltalab 2 Litre Unscrewed Brown Rectangular Graduated 24 Hour Urine Containers (White Caps, with Liner)

**Product Code:** 407003

**Brand:** Deltalab

Made of high density polyethylene, Deltalab's 2 Litre Unscrewed Brown Rectangular Graduated 24 Hour Urine Containers (White Caps, with Liner) are ideal for urine collection and other liquid samples.

They have a moulded graduation in 50ml increments up to 2 litres.

These laboratory containers are ergonomically designed for single-handed use with a lateral tap of 66mm (length).

They are supplied with an unscrewed white cap, which comes with a polexan liner to eliminate the need of a plug that will allow better handling and provide a leak-proof closure.

### Specifications

<b>Capacity (litres)</b>	2
<b>Body Colour</b>	Brown
<b>Cap Description</b>	White, with liner
<b>Presentation</b>	Unscrewed
<b>Individually Wrapped</b>	No
<b>Sterility</b>	Non-Sterile
<b>Internal Mouth (Ø mm)</b>	65
<b>Material</b>	Polyethylene
<b>Quantity</b>	43



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)