

Deltalab 2.7 Litre Unscrewed Brown Square Graduated 24 Hour Urine Containers (Non-Sterile)

Product Code: 407006

Brand: Deltalab

Made of high density polyethylene, Deltalab's 2.7 Litre Unscrewed Brown Square Graduated 24 Hour Urine Containers (Non-Sterile) are ideal for urine collection and other liquid samples.

Their ribbed caps feature a polexan liner, which eliminates the need for an insert plug and allows better handling and provides a watertight seal. These laboratory containers are graduated in 50ml increments up to 2.7 litres.

They are ergonomically designed for single-handed use and are stackable. A ring on the neck of the containers avoids dripping when pouring.



Specifications

Body Colour	Brown
Cap Description	White (with liner)
Dimensions (L x W)	130 x 126mm
Height (cap excluded)	249mm
Presentation	Unscrewed
Individually Wrapped?	No
Internal Mouth (Ø mm)	65
Material	Polyethylene
Sterility	Non-Sterile
Quantity	34

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com