

## Deltalab 2.7 Litre Unscrewed Brown Square Graduated 24 Hour Urine Containers (Non-Sterile)

**Product Code:** 407006

**Brand:** Deltalab

Made of high density polyethylene, Deltalab's 2.7 Litre Unscrewed Brown Square Graduated 24 Hour Urine Containers (Non-Sterile) are ideal for urine collection and other liquid samples.

Their ribbed caps feature a plexan liner, which eliminates the need for an insert plug and allows better handling and provides a watertight seal. These laboratory containers are graduated in 50ml increments up to 2.7 litres.

They are ergonomically designed for single-handed use and are stackable. A ring on the neck of the containers avoids dripping when pouring.

### Specifications

<b>Body Colour</b>	Brown
<b>Cap Description</b>	White (with liner)
<b>Dimensions (L x W)</b>	130 x 126mm
<b>Height (cap excluded)</b>	249mm
<b>Presentation</b>	Unscrewed
<b>Individually Wrapped?</b>	No
<b>Internal Mouth (Ø mm)</b>	65
<b>Material</b>	Polyethylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	34

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)