

Medline 5 Litre Twenty-Four Hour Urine Containers - Labelled

Product Code: 407007

Brand: Medline Scientific

Medline's 5 Litre Twenty-Four Hour Urine Containers - Labelled are rigid polyethylene jerry cans with integral handles. They feature wadded anti-leak screw caps and graduated labels with spaces for patient details.

Each urine container is CE marked in accordance with EU directive 98/79/EC.

Customer specific labelling requirements are also available.

Specifications:

| | |
|--------------------------|--------------|
| Capacity (litres) | 5 |
| Max. Width (mm) | 225 |
| Max. Depth (mm) | 180 |
| Max. Height (mm) | 280 |
| Material | Polyethylene |
| Quantity | 25 |



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com