

Deltalab 2 Litre Screwed Natural Rectangular Graduated 24 Hour Urine Containers (White Caps, with Liner) - Sterile & Individually Wrapped

Product Code: 407008.O

Brand: Deltalab

Made of high density polyethylene, Deltalab's 2 Litre Screwed Natural Rectangular Graduated 24 Hour Urine Containers (White Caps, with Liner) - Sterile & Individually Wrapped are ideal for urine collection and other liquid samples.

They have a moulded graduation in 50ml increments up to 2 litres.

These laboratory containers are ergonomically designed for single-handed use with a lateral tap of 66mm (length).

They are supplied with an unscrewed white cap, which comes with a polexan liner to eliminate the need of a plug that will allow better handling and provide a leak-proof closure.



Specifications

Capacity (L)	2
Body Colour	Natural
Cap Description	White cap, with liner
Presentation	Screwed
Individually Wrapped	Yes
Sterility	Sterile
Internal Mouth (Ø mm)	65
Material	Polyethylene
Quantity	43

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com