

Deltalab 2 Litre Screwed Natural Rectangular Graduated 24 Hour Urine Containers (Black Caps and Plug) - Individually Wrapped

Product Code: 407008OBT

Brand: Deltalab

Made of high density polyethylene, Deltalab's 2 Litre Screwed Natural Rectangular Graduated 24 Hour Urine Containers (Black Caps and Plug) - Sterile & Individually Wrapped are ideal for urine collection and other liquid samples.

They have a moulded graduation in 50ml increments up to 2 litres.

These laboratory containers are ergonomically designed for single-handed use with a lateral tap of 66mm (length).

They are supplied with a screwed white cap, which comes with a polexan liner to eliminate the need of a plug that will allow better handling and provide a leak-proof closure.

Specifications

Capacity (L)	2
Body Colour	Natural
Cap Description	Black and plug
Presentation	Screwed
Individually Wrapped	Yes
Sterility	Non-Sterile
Internal Mouth (Ø mm)	65
Material	Polyethylene
Quantity	41



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com