

Deltalab 2 Litre Non-Sterile Unscrewed Natural Graduated Cylindrical 24 Hour Urine Containers with Handles (White, with Liner and 65mm)

Product Code: 408001

Brand: Deltalab

Made of high density polyethylene, Deltalab's 2 Litre Non-Sterile Unscrewed Natural Graduated Cylindrical 24 Hour Urine Containers with Handles (White, with Liner and 65mm) are ideal for urine collection and any other liquid samples.

We recommend references with a white and ribbed cap that features a pharmaceutical grade liner, which eliminates the need for an insert plug, and allows a better handling and provides a watertight seal.

These laboratory containers are graduated in 50ml increments up to 2 litres and include a convenient handle for easy single-handed use.

Specifications

Body Colour	Natural
Cap Description	White, with liner
Base Diameter (mm)	180
Height (cap included)	290mm
Presentation	Unscrewed
Individually Wrapped?	No
Internal Mouth (Ø mm)	65
Material	Polyethylene
Sterility	Non-Sterile
Quantity	34



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com