


Dentalab 2 Litre Non-Sterile Unscrewed Natural Graduated Cylindrical 24 Hour Urine Containers with Handles (Black + Plug and 66.5mm)

Product Code: 408001OBT

Brand: Dentalab

Made of high density polyethylene, Dentalab's 2 Litre Non-Sterile Unscrewed Natural Graduated Cylindrical 24 Hour Urine Containers with Handles (Black + Plug and 66.5mm) are ideal for urine collection and any other liquid samples.

These laboratory containers are graduated in 50ml increments up to 2 litres and include a convenient handle for easy single-handed use.



Specifications

Body Colour	Natural
Cap Description	Black and plug
Base Diameter (mm)	180
Height (cap included)	290mm
Presentation	Unscrewed
Individually Wrapped?	No
Internal Mouth (Ø mm)	66.5
Material	Polyethylene
Sterility	Non-Sterile
Quantity	33

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com