


Deltalab Measuring Cylinders (25ml)

Product Code: 408025

Brand: Deltalab

Deltalab's Measuring Cylinders (25ml) have stable pentagonal bases, along with printed graduations for clear volume evaluation. Each cylinder is autoclavable up to 121°C and in compliance with ISO 6706 and BS 5404 standards.



Specifications:

Capacity (ml)	25
Graduation (ml)	5
Subdivision (ml)	0.5
Overall Diameter (mm)	19
Height (mm)	195
Material	Polypropylene
Sterility	Non-Sterile
Quantity	30

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com