


Deltalab Measuring Cylinders (100ml)

Product Code: 408100

Brand: Deltalab

Deltalab's Measuring Cylinders (100ml) have stable pentagonal bases, along with printed graduations for clear volume evaluation. Each cylinder is autoclavable up to 121°C and in compliance with ISO 6706 and BS 5404 standards.



Specifications:

| | |
|------------------------------|---------------|
| Capacity (ml) | 100 |
| Graduation (ml) | 10 |
| Subdivision (ml) | 1 |
| Overall Diameter (mm) | 31 |
| Height (mm) | 250 |
| Material | Polypropylene |
| Sterility | Non-Sterile |
| Quantity | 30 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com