

Deltalab 120ml Containers for Vacuum Systems (Individually Wrapped)

Product Code: 408726.G

Brand: Deltalab

Designed to be used with vacuum tubes, Deltalab's 120ml Containers for Vacuum Systems (Individually Wrapped) allow a clean transfer of the urine from the container to the vacuum tube.

Manufactured in transparent polypropylene and with leakproof caps made of yellow polyethylene, these laboratory containers feature an internal obturating ring. They also include a cannula with a needle inside protected by a rubber cap and a plastic suction nozzle. The caps feature a label that seals the cannula in order to prevent contact with the collection needle. The containers also feature a wide marking by the cannula.



Specifications

Max. Volume (ml, with cap)	120
Recommended Volume (ml)	60
Dimensions (mm)	57 x 73
Presentation	Individually wrapped
Sterility	Sterile
Quantity	350

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com