


Deltalab Yellow Colour Coded Inserts

Product Code: 409003AM

Brand: Deltalab

Made of polypropylene, Deltalab's Yellow Colour Coded Inserts are very convenient for colour identification in freezers or benchtops.

They fit precisely into the caps of cryovials. Each insert has a little hole in order to help remove them from the cap.



Specifications

Colour	Yellow
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com