

Deltalab Brown Polypropylene Caps for Screw Thread Microtubes (With Loop)

Product Code: 409007.M

Brand: Deltalab

Made of polypropylene, Deltalab's Brown Polypropylene Caps for Screw Thread Microtubes (With Loop) have a sealing O-ring (red) of silicone to ensure a positive leakproof seal.

For sample identification, colour coding inserts can be placed upon caps (made of polypropylene).



Specifications

Cap Dimensions (mm)	13 x 8
Description	Cap with loop
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com