


## Deltalab Transparent Polypropylene Caps for Screw Thread Microtubes (With Loop)

**Product Code:** 409007.N

**Brand:** Deltalab

Made of polypropylene, Deltalab's Transparent Polypropylene Caps for Screw Thread Microtubes (With Loop) have a sealing O-ring (red) of silicone to ensure a positive leakproof seal.

For sample identification, colour coding inserts can be placed upon caps (made of polypropylene).



### Specifications

<b>Cap Dimensions (mm)</b>	13 x 8
<b>Description</b>	Cap with loop
<b>Quantity</b>	1000

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)