

Deltalab 1.5ml Screw Thread Tamper Evident Microtubes

Product Code: 409110.2T

Brand: Deltalab

Made of autoclavable ultra clear polypropylene and with a super-fast 1/4 turn thread design, Deltalab's 1.5ml Screw Thread Tamper Evident Microtubes have ribbed caps with an internal silicone O-ring for a positive leakproof seal. Their tamper evident seal allows users to notice if the microtubes have been opened.

They are Rnase, Dnase, and Pyrogen free.


They can be used at extreme temperatures from -196°C to 121°C.

These microtubes can be used for:

- Test of fertility and DNA testing
- Packaging of diagnostic kits and reagents
- Forensic laboratories

Specifications

Volume (ml)	1.5
Microtube Dimensions	11 x 44mm
Cap Dimensions (with tamper evident ring)	15 x 9mm
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com