

Deltalab 3ml Skirted Non-Sterile Cryovials with External Threads - Graduation and Band

Product Code: 409133

Brand: Deltalab

Made of polypropylene, Deltalab's 3ml Skirted Non-Sterile Cryovials With External Threads - Graduation and Band are designed for the storage and transportation of biological material. The external thread design provides a smooth and uniform inner surface, thus reducing the risk of contamination.

The tubes can be autoclaved to 121°C and withstand temperatures down to -190°C. The closures and tubes are both manufactured of polypropylene having the same amount of expansion, which further enhances their leakproof qualities at changing temperatures.

The caps suitable for these cryovials are supplied separately.



Specifications:

Volume (ml)	3
Skirt	Yes
Dimensions (mm)	12.5 x 72
Material	Polypropylene
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com