

## Deltalab 4ml Skirted Non-Sterile Cryovials with External Threads - Graduation and and Band

**Product Code:** 409135

**Brand:** Deltalab

Made of polypropylene, Deltalab's 4ml Skirted Non-Sterile Cryovials With External Threads - Graduation and Band are designed for the storage and transportation of biological material. The external thread design provides a smooth and uniform inner surface, thus reducing the risk of contamination.

The tubes can be autoclaved to 121°C and withstand temperatures down to -190°C. The closures and tubes are both manufactured of polypropylene having the same amount of expansion, which further enhances their leakproof qualities at changing temperatures.

The caps suitable for these cryovials are supplied separately.

### Specifications:

<b>Volume (ml)</b>	4
<b>Skirt</b>	Yes
<b>Dimensions (mm)</b>	12.5 x 76
<b>Material</b>	Polypropylene
<b>Quantity</b>	1000

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)