

Deltalab 0.5ml Microtubes (Natural Graduated)

Product Code: 4092.IN

Brand: Deltalab

Made of ultra clear polypropylene and autoclavable, Deltalab's 0.5ml Microtubes (Natural Graduated) have a frosted writing surface on the side and on closure for better sample identification. The caps are easy to handle and can be opened and closed with one hand. The lids are easily pierceable. These microtubes have moulded-in graduations in 100µl increments from 0.1 to 0.6ml. They can withstand centrifugation up to 17,000xg and be used from -80°C to 121°C.

They are RNAse, DNAse, and PCR inhibitors free.



Specifications

Dimensions (mm)	30 x 8
Description	Natural graduated
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Features: