

Deltalab 1.5ml Graduated Microtubes

Product Code: 4092.2N

Brand: Deltalab

Made of ultra clear polypropylene and autoclavable, Deltalab's 1.5ml Graduated Microtubes have a frosted writing surface on the side and on closure for better sample identification. The caps are easy to handle and can be opened and closed with one hand. The lids are easily pierceable and be used from -170°C to 121°C.

The caps are specially designed to reduce popping when centrifugation, freezing, boiling, storing or shipping is required.

When the objective is to avoid liquid remaining on the wall of the tube, just push the cap slightly and a little centrifugation time is required.

When centrifugation is required, push the cap firmly to achieve complete closure.

The tubes can withstand centrifugation up to 20,000xg.



Specifications

Dimensions (mm)	39 x 11
Colour	Natural
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com