

## Deltalab 1.5ml DNase & RNase Microtubes


**Product Code:** 4092.3N

**Brand:** Deltalab

Made of ultra clear polypropylene, Deltalab's 1.5ml DNase & RNase Microtubes are autoclavable and feature attached flat caps. They are graduated from 100µl.

The tubes can be easily opened and closed with one hand. They are ideal for the storage and freezing of any kind of biological material (cells, sperm, bacteria, etc.).

They can be used at temperatures from -90°C to 121°C and withstand centrifugation up to 20,000xg.



### Specifications

<b>Dimensions (Ø x h; mm)</b>	10.7 x 39.1
<b>Description</b>	Natural colour
<b>Quantity</b>	500

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)