

Detalab 2ml Microtubes with Attached Caps, Natural Colour & Clear Siliconised

Product Code: 4092.6NS

Brand: Detalab

Made of ultra clear polypropylene, Detalab's 2ml Microtubes with Attached Caps, Natural Colour & Clear Siliconised feature a frosted writing surface on the side and on closure for better sample identification. They are autoclavable and certified RNase and DNase free.


The tubes can be opened and closed with one hand. The lids are easily pierceable. They have moulded-in graduations in 100µl increments up to 2ml.

They can be used at temperatures from -90°C to 121°C and can withstand centrifugation up to 17,000xg.

They also have a low adhesion surface for special applications such as protein work.

Specifications

Dimensions (mm)	9 x 40
Description	Natural colour, clear siliconised
Quantity	250



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com