

## Deltalab 2ml Microtubes with Attached Caps, Natural Colour & Clear Siliconised

**Product Code:** 4092.6NS

**Brand:** Deltalab

Made of ultra clear polypropylene, Deltalab's 2ml Microtubes with Attached Caps, Natural Colour & Clear Siliconised feature a frosted writing surface on the side and on closure for better sample identification. They are autoclavable and certified RNase and DNase free.


The tubes can be opened and closed with one hand. The lids are easily pierceable. They have moulded-in graduations in 100µl increments up to 2ml.

They can be used at temperatures from -90°C to 121°C and can withstand centrifugation up to 17,000xg.

They also have a low adhesion surface for special applications such as protein work.

### Specifications

<b>Dimensions (mm)</b>	9 x 40
<b>Description</b>	Natural colour, clear siliconised
<b>Quantity</b>	250



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)